

Hauf Lab Tools – Quick Start Guide

Installation

Ensure that your ImageJ/FIJI is of the latest version

Step 1. In Fiji, go to **Help > Update...**

Step 2. In the ImageJ Updater, click **Manage update sites**

Step 3. In the Manage update sites window, click **Add update site**
A new update site row should appear at the bottom of the list


Step 4. Under the Name column, change **New** to **Hauf Lab**
Under the URL column, add **https://sites.imagej.net/HaufLab**

Step 5. Make sure that the Hauf Lab update site row is checked in the far-left column.
You can close the Manage update sites window with the **Close** button.

Step 6. In the ImageJ Updater, click **Apply changes**

Step 7. Restart Fiji.

Usage

Step 1. In the Fiji toolbar, click the double-arrow  button.

Step 2. From the menu, click **Hauf_Lab_Tools**
This initializes the Hauf Lab Tools menu – you will need to do this each time you start Fiji and wish to use Hauf Lab Tools. This also updates the macro list of Hauf Lab Tools to match those that are locally stored.

Step 3. A new icon should appear in the Fiji toolbar, the **Hauf Lab Tools Menu** 

Step 4. From the Hauf Lab Tools Menu, select your desired macro to run

Source Code

To easily see the source code of any of the macros from Hauf Lab Tools, currently installed on your system, enable the built in 'Additional Actions Prompts'.

In the **Hauf Lab Tools Menu** , go to **Options**

Enable the **Additional Actions Prompt** by checking the checkbox.

Click **OK** to apply changes.

Now whenever you select a macro from the Hauf Lab Tools Menu, an additional prompt will appear, allowing you to either run your selected macro, or open the locally-stored source code.